

A&D MEDICAL THERMOMETER

FROM:

AD
A & D MEDICAL



A lightweight and compact thermometer that provides quick and accurate temperature measurements at the push of a button.

The A&D Instant Read Ear Thermometer has a unique and attractive design that provides a fast, accurate and easy way to take a temperature through the ear. A temperature reading is provided in just one second! The Thermometer does not require a protective lens cover—its waterproof tip can be easily cleaned with water or rubbing alcohol.

- An accurate reading in just one second (scans the ear 512 times per second.)
- Waterproof lens tip that can be easily cleaned with water or rubbing alcohol.
- Stores the last reading in memory.
- Automatic power-off after 1 minute.
- Readings in Fahrenheit or Celsius.
- Includes protective carrying case for safe storage.
- Trilingual instruction manual with illustrations and quick start card in English, French and Spanish.

A&D Medical Pty Ltd

Ph: 1800 01 112

Email: 1800 011 112

Web: www.andmedical.com.au