



UPPER ARM BLOOD PRESSURE MONITOR

Model: UA-611



KEY FEATURES:

-  One Button Operation
-  Irregular Heartbeat Detection
-  Blood Pressure Rating Indicator
-  30 Reading Memory
-  5 Year Warranty
-  Cuff
22-32 cm

A classic blood pressure monitor with the benefits of speed and accuracy at a touch of a button.

Trusted Accuracy

- Clinically validated for accuracy
- Indicates an irregular heartbeat (IHB) during a measurement
- Pressure Rating Indicator classifies blood pressure in accordance with medically approved standards

User Friendly

- Easy one button operation for quick and accurate measurement
- Stores 60 measurements in memory with recall
- Displays average of all measurements
- Gentle Inflation enables a more comfortable measurement

The A&D Promise

- 5 year warranty on the monitor / 2 years on the cuff
- Cuffs are not made with natural rubber latex
- Battery life indicator
- Download the A&D Heart Track App for Free to manually log and electronically track measurements (available on iOS and Android)

Technical Specifications

Measurement Method	Oscillometric
Memory	30 measurements
Cuff Size	22 - 32cm adult cuff size
Measurement Range	Pressure: 0-299 mmHg Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg Pulse: 40-180 beats/min
Measurement Accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA)
Battery Life	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
Monitor Weight	230 g without batteries
Operating Conditions	+10°C to +40 °C / 15 %RH to 85 %RH / 800 hPa to 1060 hPa
Storage / Transport Conditions	-20 °C to +60 °C / 10 %RH to 95 %RH 700hPa to 1060hPa

Packaging and shipping specifications

Dimensions	Width	Height	Diameter	SKU
Monitor Dimensions	96 mm	68 mm	130 mm	09010016

What is In The Box

- 1 Monitor
- 1 Blood pressure cuff (22cm to 32cm)
- 1 Instruction Manual

