



PRECISION PERSONAL HEALTH SCALE

Model: UC-321







A&D's outstanding technology for weighing instruments introduces the low profile UC-321 for body weight measurement at home. A high resolution of 50g provides ultra precise weighing results in the comfort of your own home.

General Specifications

- Maximum capacity of 150kg
- 50g resolution provides ultra precise weighing
- Built-in memory function can store up to 31 readings in its memory that can be recalled with a press of a button
- Target weight setting with BMI value calculations
- Ultra low profile so its easy to step onto the scale
- LCD character height: 25mm
- Battery operated with 4 x AA cell batteries
- Battery life is approximately 2,000 measurements

UC-321 SCALE SPECIFICATIONS	
Capacity	150kg
Resolution	50g

^{*} Specifications are subject to change without notice

