

A&D Medical ... *Clearly a Better Value*

## PRECISION HEALTH SCALE

LARGE PLATFORM SIZE FOR BETTER PATIENT STABILITY - 360mm x 415mm



- High Capacity - 250kg (550 lb)
- Low Profile
- Step-on technology
- Large Display
- Includes Scale tray
- 200 measurement memory

The UC-355PBT-Ci Precision Health Scale has a 250kg capacity, supporting a large range of users. This scale can display your weight in kilograms or pounds, allowing you to select what you are comfortable with. Additionally, the scale has a wider platform to support those with a wide base. It is also equipped with audible buzzers and visual indicators to let you know when to get on and off. The profile is less than 35mm tall even when the carpet tray (included) is utilised. You don't have to push any buttons to turn on this scale, with A&D's step-on technology. Just step on and let the scale guide you with on-screen prompts.

# PRECISION HEALTH SCALE

PART NUMBER	DESCRIPTION
09030901	250kg Precision Health Scale



## Technical Specifications

Measuring method	Electric resistance
Digital Display	Liquid crystal display, 35 mm
Maximum capacity	250 kg / 550 lbs
Resolution	0.1 kg between 0 - 100 kg / 0.2 lb between 0 - 220 lb 0.2 kg between 100 - 250 kg / 0.5 lb between 220 - 550 lb
Power Source	Four type AA batteries (1.5 volt)
Dimensions	360 mm (d) x 415 mm (w) x 39 mm (h)
Weight	Approximately 4.6 kg, including the batteries
Accessories	Instruction manual, batteries (AA x 4), scale tray
Battery Life	Approximately 1 year with two measurements per day
Setting Items (Units)	kg or lb

Ask your A&D Medical Sales Representative about other A&D Medical products.

- Blood Pressure Monitors
- Digital Thermometers
- Personal Scales
- Activity Monitors

For questions about A&D products, call us  
1800 011 112