



AD-6121A

Precision Scale

The AD-6121A precision scale was designed for professional use with slim and compact functions in mind. A low profile of only 57 mm makes it easily accessible for patients and its 100 g resolution can display small differentiations in weight. This precision scale weighs only 3 kg so it can be easily carried to a suitable location for patient weighing and its compact size means that it can be easily stored under a patient's bed.



Features

- EC Type approved Class III
- Standard uni-directional RS-232C
- 40-measurement memory
- 100 g resolution and up to 150 kg capacity
- Very low profile (57 mm height) with rugged construction
- Light as 3 kg for great portability
- Built-in handle for easy transportation
- Back-lit display for easier viewing in dim conditions
- Tare function for precise weighing
- Measurement values can be rotated 180 degrees for ease of use during group medical checkups

Specifications

Capacity	150 kg
Resolution	100 g
Sensor	Load cell
Display	LCD, character height: 25mm
Power	4 type AA batteries (R6P)
Battery life	Approx. 2,000 measurements
Weight	Approx. 2.7 kg
Dimensions	320 x 314 x 35 mm (WxDxH)
Auto power off	2 minutes
Operating temperature	5 to 35°C (50 to 95°F)
Communication	EIA RS-232C (uni-directional)

AND
A&D Medical
<http://www.aandd.jp>