



# HEIGHT MEASUREMENT SCALE

## Model: AMP-250

### General Specifications

- Maximum capacity of 250kg.
- 100g resolution provides precise weighing.
- BMI function to determine underweight/ normal weight / surplus weight.
- Bright LCD display with white colour LED backlight.
- Easy to use height measuring stick (88-205cm).
- Non slip base with rubber feet to provide secure footing.
- Lock and Hold functions are available as standard.
- Easy to transport with two integrated transport wheels.
- Platform size is 370mm x 370mm.
- Powered by AC adaptor with a rechargeable battery.
- Operating time up to 40 hours without backlight, and the charging time is approximately 12 hours.



#### AMP-250 SCALE SPECIFICATIONS

Capacity	250kg
Resolution	100g

\* Specifications are subject to change without notice