

ERGONOMICALLY DESIGNED CHAIR SCALE

Model: HVL-CS

General Specifications

- Simple operation with 2 button keypad
- Easy to clean chair shell to ensure maximum hygiene
- Long battery life (up to 600hrs) on alkaline cells
- Highly portable with locking castors for patient safety
- High visibility display with easy viewing from multiple angles and auto power off after 3 minutes of inactivity
- Ergonomically designed so minimal effort is required by patient and staff
- Push handles make it convenient for staff when moving the scale
- Australian made, quality manufacturing (ISO-9001)



HVL-CS SCALE SPECIFICATIONS

Capacity	Up to 150kg
Minimum resolution	20g
Weight display	Liquid Crystal Display (LCD)
Power	6 x D cell batteries or AC adaptor (optional)
Chair shell material	White high strength moulded plastic
Base material	Grey die cast aluminium spider
Wheels dimension	5 x 125mm diameter
Safety	All wheels are fitted with locking brakes
Loadcell	Allows for up to 200% overload
Overall dimensions (mm)	750 (w) x 630 (d) x 880 (h)
Packed dimensions (mm)	930 (w) x 670 (d) x 950 (h)

* Specifications are subject to change without notice