

HIGH CAPACITY PERSONAL SCALE


Model: PW-200-FG

General Specifications

- High capacity platform scale (up to 200kg)
- Excellent resolution accuracy (20g)
- Large non-slip weigh platform
- Powered by AC adaptor or 4 x C cell batteries
- Easy Portability
- Clear LCD display, and easy to use keyboard
- RS-232C interface option available



Dimensions (mm)



PW-200-FG SCALE SPECIFICATIONS

Capacity	200kg
Resolution	20g

* Specifications are subject to change without notice