

# **WCS**

# **Wheelchair Scale**



- **300 kg Capacity**
- **Lightweight - easy to transport manually**
- **Simple Operation**
- **Reduces Patient Disturbance**
- **Ergonomically Designed**

**A&D**

**A&D Australasia**

*"Excellence in Measurement"*

[www.andaustralasia.com.au](http://www.andaustralasia.com.au)





## Features

### **Lightweight Drop Centre Design**

Easy access to drop platform surface via entry ramp

### **Anodised Finish**

Visually appealing for your location

### **Dual Power Supply**

Unit can be operated on AC or battery power

### **300 kg Capacity - 100 gram Accuracy**


Suitable for all wheelchairs

### **Separate Indicator Stand**

Allows maximum flexibility in set up

### **Large Weights Area**

Suitable for both motorised and non-motorised wheelchairs



**A&D Australasia**  
*"Excellence in Measurement"*  
 ACN 007 556 809  
 Head Office  
 32 Dew Street  
 THEBARTON South Australia 5031  
 Telephone (08) 8301 8100  
 Facsimile (08) 8352 7409

Victorian Office  
 39 Bakehouse Road  
 KENSINGTON Victoria 3031  
 Telephone (03) 9372 1522  
 Facsimile (03) 9372 1193

New South Wales Office  
 Unit 4, 14 Abbott Road  
 SEVEN HILLS New South Wales 2147  
 Telephone (02) 9674 5466  
 Facsimile (02) 9674 2544

Email: [sales@andaustralasia.com.au](mailto:sales@andaustralasia.com.au)  
 Internet page: [www.andaustralasia.com.au](http://www.andaustralasia.com.au)

Your Dealer