

PREMIUM HIGH CAPACITY PERSONAL SCALE

Model: UC-356BLE

A&D's outstanding technology for weighing instruments introduces the UC-356BLE. A wide base personal weight scale with high-capacity up to 250kg.



General Specifications

- High capacity of 250kg (100kg x 100g/100-250kg x 200g)
- Stable footprint, Dimensions (40 x 415 x 361 mm)
- 4 precision sensors to ensure accurate readings
- Motion technology for accurate measurements even with slight movement
- Measurements are recorded in pounds or kilograms
- User Friendly, step on and measure your weight.
- LCD display for easy reading
- Stores 63 measurements
- Auto off to save battery life
- 2 year warranty

Technical Specifications

Maximum capacity	250kg
Internal Memory	63 measurements
Minimum display	0 to 100kg with an increment of 0.1 kg, and 100 to 250kg with an increment of 0.2 kg
Wireless communication	Bluetooth Ver.4.0LE
Units	kg or lb
Power source	Four AA size alkaline batteries
Battery life	Approx. 1 year (two measurements per day when alkaline batteries are used)
Dimensions (HxWxD)	40.6 x 414.8 x 360.7 mm
Weight	3.5kg
Data Outputs via Integration	Time, date, weight value, battery status, set memory

* Specifications are subject to change without notice