

DIGITAL BLOOD PRESSURE MONITOR UB-511

Clinically
Validated

IHB



Fast and Ultra Compact Wrist Monitor

- Ultra Compact and Light Weight Design
- Fast Measurement along with Inflation
- Irregular Heart Beat (IHB) Indicator
- Blood Pressure Classification Indicator
- 30 Memory + Average Reading

CE 0123

A&D
A&D Medical

Technology You Can Trust

<http://www.aandd.jp>

DIGITAL BLOOD PRESSURE MONITOR UB-511

Fast and Ultra Compact Wrist Monitor

Ultra Compact and Light Weight Design

Our state-of-the-art technology brings real portability to the wrist monitor. It is 40% smaller and 30% lighter than our former model. Such a miniature monitor promises portable measurement whenever user needs to monitor his/her blood pressure.

Fast Measurement Along with Inflation

The UB-511/512 series completes blood pressure measurement while inflating the wrist cuff. This gives user a quick measurement anywhere he/she needs to monitor his blood pressure.

Irregular Heart Beat (IHB) Indicator

This lights up if IHB is detected during the measurement. Helps users to find an irregular heartbeat which may be associated with cardiovascular disease. (Note: This is for warning purposes only. A&D recommends contacting a physician if the user sees this indicator frequently)

Blood Pressure Classification Indicator


World Health Organization (WHO) classifies blood pressure ranges into 6 stages. This monitor indicates your position (Optimal to Severe Hypertension) automatically after each measurement is complete.

30 Memory + Average Reading

The UB-511 has large memory of 30 readings for monitoring the trend of measurements. An average of all the stored data is shown at each start up, giving a condensed reading history at a glance.

WHO blood pressure classification

Standards for assessment of high blood pressure, without regard to age, have been established by the World Health Organization (WHO), as shown in the chart on the right.



A&D
A&D Medical

Technology You Can Trust

<http://www.aandd.jp>

*UB511CA-ADCC-03-PR3-143**
REF*CARD-UB511